# HISTO TRAY RESERVOIR

## HLA Klasse I

**IN VITRO DIAGNOSTIKA**

<b>HISTO TRAY AB 72 RESERVOIR</b>		<b>6283</b>	<b>HISTO TRAY B27 (20) RESERVOIR</b>		<b>6291</b>
<b>HISTO TRAY ABC 72 RESERVOIR</b>		<b>6288</b>	<b>HISTO TRAY B27 (50) RESERVOIR</b>		<b>6290</b>
<b>HISTO TRAY ABC 120 RESERVOIR</b>		<b>6286</b>			
<b>HISTO TRAY ABC 144 RESERVOIR</b>		<b>6287</b>			

**Produktbeschreibung**

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**Description of product**

**HISTO TRAY RESERVOIR** are storage trays with anti-HLA sera and controls for tissue typing of HLA class I antigens. Worksheets for evaluation and the listing of test results are enclosed.

**Test principle**

Anti-HLA sera react with the correspondent membrane-bound antigens on human lymphocytes. The addition of rabbit complement results in a structural change of the cell membrane which leads to a penetration of an indicator dye. Stained lymphocytes = positive reaction. In case of missing antigen-antibody reaction, the cell membrane is intact. No penetration of indicator dye takes place and the cells remain unstained = negative reaction.

**Additional materials required**

Aqua dest., paraffin perliquidum

Microtest plates with 60 or 72 wells (e.g. Greiner), automatic dropper for 1 – 5 µl, microliter syringe for 1 – 5 µl

Optional: Vacuum packing machine, vacuum bags

**NIH technique**

Cell culture medium (e.g. RPMI 1640), cell separation medium (e.g. HISTOPREP), centrifuge tubes (12 ml), pasteur pipettes, centrifuge, Neubauer count chamber or cell counter

Rabbit complement (10 x 1 ml **REF** 7018, 10 x 5 ml **REF** 7023)

Eosin solution (5% aqueous), Formaldehyde solution (37%, pH 7.2)

Cover glasses, inverse phase contrast microscope

**Immuno Beads method**

Immuno Beads (details see manufacturer)

Rabbit complement

Acridinorange/Ethidiumbromide, quenching-solution, EDTA (8% aqueous)

Fluorescence microscope

**Test procedure – Dropping of microtest plates**

1. Thaw storage tray with anti-HLA sera and controls.
2. Label microtest plates with product name, lot and expiry date (water proof labelling).
3. Flush microliter syringe or automatic dropper at least 3 x with fresh aqua dest.. If using an automatic dropper, note the instruction manual of the manufacturer.
4. Drop 5 µl paraffin perliquidum in each well of the microtest plate.
5. Drop 1 µl anti-HLA serum or control from the storage tray into the respective well.  
Pay attention that the red coloured serum is located under the paraffin perliquidum at the bottom of the well.
6. Store the dropped microtest plates at ≤ -20°C up to use.  
It is recommended to weld the dropped plates in vacuum bags under vacuum for storage.

**Test procedure - Isolation of lymphocytes from e.g. heparinized blood**

1. In order to increase the cell yield, dilute 4 ml of heparinized blood (50 I.U./ml) with 4 ml of cell culture medium (e.g. RPMI 1640).
2. Pipet 4 - 5 ml of cell separation medium, e.g. **HISTOPREP** into a centrifuge tube (12 ml).
3. Carefully add app. 6 ml of diluted blood with a Pasteur pipette to the gradient alongside the inner edge of the tube.
4. Centrifuge 15 minutes at 1200 x g and a temperature of 18...22°C (room temperature). Centrifuge without braking.
5. Take off the lymphocyte ring (interphase) using a Pasteur pipette and pipet it into a new centrifuge tube.
6. For lymphocyte washing, fill it up with cell culture medium, e.g. RPMI 1640, and centrifuge for 10 minutes at 550 x g. Discard the supernatant, resuspend the sediment and fill it up with cell culture medium, e.g. RPMI 1640.
7. Centrifuge for 10 minutes at 230 x g, discard the supernatant, resuspend the bottom sediment and fill it up with cell culture medium, e.g. RPMI 1640.
8. Centrifuge for 10 minutes at 110 x g and discard the supernatant.
9. Resuspend the sediment in cell culture medium, e.g. RPMI 1640, and adjust to a final concentration of 2000 - 3000 lymphocytes per µl (Neubauer count chamber or cell counter).

**Test procedure – NIH technique**

1. Bring the microtest plates with predropped anti-HLA sera to a temperature of 18...22°C (room temperature).
2. Place 1 µl lymphocyte suspension (2.000 - 3.000 cells) into each predropped well.  
In order to guarantee sufficient antigen-antibody reaction it is necessary that antiserum and cells touch each other.
3. Incubate at a temperature of 18...22°C (room temperature) for 30 minutes.
4. Add 5 - 6 µl rabbit complement.
5. Incubate at a temperature of 18...22°C (room temperature) for 60 minutes.
6. Add 3 - 4 µl Eosin solution (5% aqueous) (soft touch method) and incubate for 5 - 10 minutes.
7. Fix with 5 - 6 µl Formaldehyde solution (37%, pH 7.2) (soft touch method). Allow sedimentation of cells at least 60 minutes.
8. Cover the tray with a cover glass shortly before reading under an inverse phase contrast microscope.

**Test procedure - Isolation of T-lymphocytes from e.g. heparinized blood**

Isolation of the T-lymphocytes using the Immuno Beads method as well as staining and fixation reagents according to manufacturers instructions.

**Test procedure - Immuno Beads (IMB) method**

1. Bring the microtest plates with predropped anti-HLA sera to a temperature of 18...22°C (room temperature).
2. Place 1 µl IMB-T-lymphocyte suspension (app. 1.000 cells) into each predropped well.  
In order to guarantee sufficient antigen-antibody reaction it is necessary that antiserum and cells touch each other.
3. Incubate at a temperature of 18...22°C (room temperature) for 30 minutes.
4. Add 5 µl rabbit complement Acridinorange/Ethidiumbromide (AO/EB) (1.000 µl rabbit complement + 20 µl AO/EB).
5. Incubate for 60 minutes at a temperature of 18...22°C (room temperature) in darkness.
6. Add 5 µl EDTA/quenching-solution (2.000 µl quenching solution + 1.000 µl EDTA 8% aqueous).
7. Read under a fluorescence microscope.

## Evaluation of results

The amount of lysed lymphocytes compared with the total amount of lymphocytes is quoted as a score value in each well.

% lysed cells		Score	Evaluation
0	- 19%	= Score 1	negative
20	- 39%	= Score 2	doubtful negative
40	- 59%	= Score 4	weak positive
60	- 79%	= Score 6	positive
80	- 100%	= Score 8	strong positive
		= Score 0	no evaluation possible

## Troubleshooting

### Causes of false negative or weak reactions

- Erythrocyte contamination can make microscopic evaluation difficult
- Platelet contamination
- The amount of lymphocytes is too high
- The amount of anti-HLA serum is too low
- Yellow colour of the HLA antisera
- Trays have been thawed and refrozen
- Reconstituted complement kept too long at room temperature before use
- Residual complement was frozen and thawed again
- Incubation time were too short
- Incubation temperature were too low

### Causes of false positive reactions

- Cross reactions
- The amount of anti-HLA serum is too high
- Incubation time were too long
- Incubation temperature were too high
- Prior damage of lymphocytes (negative control is positive = „background“)
- Failure to add fixative

## Rabbit Complement

Reconstitute lyophilized complement with 1 ml or 5 ml aqua dest. (according to package sizes). The reconstitution takes 10 - 15 minutes. Reconstituted complement must be stored cool (2...8°C) and used within 3 - 4 hours.

**DO NOT FREEZE** dissolved rabbit complement!

## Performance characteristics

Please note the listing of test results HISTO TRAY RESERVOIR in order to receive data for diagnostic sensitivity and specificity (R-Value).

## Literature

Bodmer, J. et al., 1997. Tissue Antigens 49:297-321

## Warnings and Precautions

**HISTO TRAY RESERVOIR** are designed for in vitro diagnostic use only and should be applied by properly trained personnel, experienced in histocompatibility testing. Transplantation guidelines as well as EFI standard should be followed, in the particular case of doubtful typing results. Human source material used to produce these reagents has been tested and found negative for HBsAg and HIV and HCV antibodies. Nevertheless all used biological material like blood, sera and control sera should be handled as potentially infectious, because no test method can guarantee that material derived from biological sources are free from infectious agents. When handling biological material appropriate safety precautions are recommended (Do not pipet by mouth; wear disposable gloves while handling biological material and performing the test; disinfect hands when finished the test). Biological material should be inactivated before disposal (e.g. in an autoclave). Disposables should be autoclaved or incinerated after use. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated areas swabbed with a suitable standard disinfectant or 70% alcohol. Material used to clean spills, including gloves, should be inactivated before disposal (e.g. in an autoclave).

Anti-HLA sera contain as preservative < 0.1% NaN<sub>3</sub>. A concentration of < 0.1% NaN<sub>3</sub> is not considered to be a harmful concentration. Nevertheless avoid contact with the skin and mucous membranes. The copper and lead used in some plumbing systems can react with azides to form explosive salts. The quantities of azide used in this reagent are small; nevertheless when disposing of azide-containing materials, they should be flushed away with a large volume of water.

Disposal of all specimen and test materials should be in accordance with state and local law.

For quenching solution, Formaldehyde solution and Acridinorange /Ethidiumbromide (AO/EB) please note the warnings and precautions of the manufacturers.

A yellow colouration of anti-HLA sera which still remains after thawing, may indicate a change of the pH value. Those plates should **not** be used for the test.

Material Safety Data Sheets available on request.

Do not use **HISTO TRAY RESERVOIR** beyond the indicated expiration date on the label.

**Preservative:** < 0.1% NaN<sub>3</sub>

**Storage:** ≤ - 20°C

**Shelf life:** until the expiration date indicated on the labels

**Package:** according to information indicated on the labels

## Explanation of symbols used on Labelling

<b>IVD</b>	For in vitro diagnostic use		Storage temperature
<b>HLA TYPING</b>	Intended purpose: HLA typing		Use by
<b>LOT</b>	Batch code		Consult Instructions for use
<b>REF</b>	Catalogue number		

**Version: 1 / 2015 | Issue: 2015-02**

**BAG Health Care GmbH**

Amtsgerichtsstraße 1-5  
35423 Lich / Germany

Tel.: +49 (0) 6404 / 925 - 0  
Fax: +49 (0) 6404 / 925 - 250

www.bag-healthcare.com

info@bag-healthcare.com

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Tel.: +49 (0) 6404 / 925 - 450

Fax: +49 (0) 6404 / 925 - 460

verkauf@bag-healthcare.com

**Customer Service:**

Tel.: +49 (0) 6404 / 925 - 125

Fax: +49 (0) 6404 / 925 - 421

service@bag-healthcare.com