## HISTO TRAY HLA Klasse I

|                               |            |             |                                  |            |             |
|-------------------------------|------------|-------------|----------------------------------|------------|-------------|
| <b>HISTO TRAY AB 72 (10)</b>  | <b>REF</b> | <b>7027</b> | <b>HISTO TRAY B27 (10)</b>       | <b>REF</b> | <b>7006</b> |
| <b>HISTO TRAY AB 120 (5)</b>  | <b>REF</b> | <b>7042</b> | <b>HISTO TRAY B27 (20)</b>       | <b>REF</b> | <b>7008</b> |
| <b>HISTO TRAY AB 144 (5)</b>  | <b>REF</b> | <b>7043</b> | <b>HISTO TRAY B27 (50)</b>       | <b>REF</b> | <b>7009</b> |
| <b>HISTO TRAY ABC 72 (10)</b> | <b>REF</b> | <b>7022</b> | <b>HISTO TRAY B27 forte (10)</b> | <b>REF</b> | <b>7004</b> |
| <b>HISTO TRAY ABC 120 (5)</b> | <b>REF</b> | <b>7013</b> |                                  |            |             |
| <b>HISTO TRAY ABC 144 (5)</b> | <b>REF</b> | <b>7035</b> |                                  |            |             |

### Produktbeschreibung

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### Description of Product

**HISTO TRAY kits** contain microplates with predropped anti-HLA sera and controls. Rabbit complement, worksheets for evaluation and the listing of test results are enclosed. The intended use of **HISTO TRAY AB-** and **-ABC kits** is the serological tissue typing of HLA class I antigens. The **HISTO TRAY B27 kits** are used for the determination of HLA-B27. Based on significant associations between HLA-B27 and a group of diseases summarized as seronegative arthritis (Bechterew's disease, Reiter's disease, reactive arthritis) the evidence of this antigen will be used to support the diagnosis of these diseases.

### Test Principle

HLA-antisera react with the corresponding membrane-bound antigens on human lymphocytes. The addition of rabbit complement results in a structural change of the cell membrane which leads to a penetration of an indicator dye. Stained lymphocytes = positive reaction. In case of missing antigen-antibody reaction, the cell membrane is intact. No penetration of indicator dye takes place and the cells remain unstained = negative reaction.

### Test Procedure - Isolation of Lymphocytes from e.g. Heparinized Blood

- In order to increase the cell yield, dilute 4 ml of heparinized blood (50 I.U./ml) with 4 ml of cell culture medium (e.g. RPMI 1640).
- Pipet 4 - 5 ml of cell separation medium, e.g. HISTOPREP, into a centrifuge tube (12 ml).
- Carefully add app. 6 ml of diluted blood with a Pasteur pipette to the gradient alongside the inner edge of the tube.
- Centrifuge 15 minutes at 1200 x g and a temperature of 18...22°C (room temperature). Phase out centrifugation without using brake function.
- Take off the lymphocyte ring (interphase) using a Pasteur pipette and pipet it into a new centrifuge tube.
- For lymphocyte washing, fill it up with cell culture medium, e.g. RPMI 1640, and centrifuge for 10 minutes at 550 x g. Discard the supernatant, resuspend the sediment and fill it up with cell culture medium, e.g. RPMI 1640.
- Centrifuge for 10 minutes at 230 x g, discard the supernatant, resuspend the bottom sediment and fill it up with cell culture medium, e.g. RPMI 1640.
- Centrifuge for 10 minutes at 110 x g and discard the supernatant.
- Resuspend the sediment in cell culture medium, e.g. RPMI 1640, and adjust to a final concentration of 2000 - 3000 lymphocytes per µl (Neubauer count chamber or cell counter).

### Test Procedure – NIH Technique

- Bring the HISTO TRAY plates to a temperature of 18...22°C (room temperature).
- Place 1 µl lymphocyte suspension (2.000 - 3.000 cells) into each predropped well.  
In order to guarantee sufficient antigen-antibody reaction it is necessary that antiserum and cells touch each other.
- Incubate at a temperature of 18...22°C (room temperature) for 30 minutes.
- Add 5 - 6 µl rabbit complement.
- Incubate at a temperature of 18...22°C (room temperature) for 60 minutes.
- Add 3 - 4 µl Eosin solution (5% aqueous) (soft touch method) and incubate for 5 - 10 minutes.
- Fix with 5 - 6 µl Formaldehyde solution (37%, pH 7.2) (soft touch method). Allow sedimentation of cells at least 60 minutes.
- Cover the tray with a cover glass shortly before reading under an inverse phase contrast microscope.

### Test Procedure - Isolation of T-lymphocytes with the Immunobeads Technique

Please refer to the manufacturer's instructions

- when using the Immunobeads (IMB) technique for isolation of the T-lymphocytes
- regarding the reagents needed for staining and fixation

### Test procedure - Immunobeads (IMB) Technique

- Bring the HISTO TRAY plates to a temperature of 18...22°C (room temperature).
- Place 1 µl IMB-T-lymphocyte suspension (app.1.000 cells) into each predropped well.  
In order to guarantee sufficient antigen-antibody reaction it is necessary that antiserum and cells touch each other.
- Incubate at a temperature of 18...22°C (room temperature) for 30 minutes.
- Add 5 µl rabbit complement Acridinorange/Ethidiumbromide (AO/EB) (1.000 µl rabbit complement + 20 µl AO/EB).
- Incubate for 60 minutes at a temperature of 18...22°C (room temperature) in darkness.
- Add 5 µl EDTA-/quenching-solution (2.000 µl quenching solution + 1.000 µl EDTA 8% aqueous).
- Read HISTO TRAY plates under a fluorescence microscope.

### Evaluation of Results

The amount of lysed lymphocytes compared with the total amount of lymphocytes is quoted as a score value in each well.

| % lysed cells | = | Score   | Evaluation              |
|---------------|---|---------|-------------------------|
| 0 - 19%       | = | Score 1 | negative                |
| 20 - 39%      | = | Score 2 | doubtful negative       |
| 40 - 59%      | = | Score 4 | weak positive           |
| 60 - 79%      | = | Score 6 | positive                |
| 80 - 100%     | = | Score 8 | strong positive         |
|               | = | Score 0 | evaluation not possible |

### Troubleshooting

#### **Causes of false negative or weak reactions**

- Erythrocyte contamination can make microscopic evaluation difficult
- Platelet contamination
- The amount of lymphocytes is too high

- Yellow colour of the HLA antisera
- Trays have been thawed and refrozen
- Reconstituted complement kept too long at room temperature before use
- Residual complement was frozen and thawed again.
- Incubation times were too short
- Incubation temperatures were too low

#### Causes of false positive reactions

- Cross reactions
- Incubation times were too long
- Incubation temperatures were too high
- Prior damage of lymphocytes (negative control is positive = „background“)
- Failure to add fixative

#### Rabbit Complement

Reconstitute lyophilized complement with 1 ml aqua dest.. The reconstitution takes 10 - 15 minutes. Reconstituted complement must be stored cool (2...8°C) and used within 3 - 4 hours.

**DO NOT FREEZE** dissolved rabbit complement!

#### Performance Characteristics

Please refer to the listing of HISTO TRAY test results to receive data for diagnostic sensitivity and specificity (R-Value).

#### Literature

Bodmer, J. et al., 1997. Tissue Antigens 49:297-321

#### Warnings and Precautions

**HISTO TRAY plates and rabbit complement** are designed for in vitro diagnostic use only and should be applied by properly trained personnel, experienced in histocompatibility testing. Transplantation guidelines as well as EFI standard should be followed, especially in the particular case of doubtful typing results. Human source material used to produce these reagents has been tested and found negative for HBsAg and HIV and HCV antibodies. Nevertheless all used biological material like blood, sera and control sera should be handled as potentially infectious, because no test method can guarantee that material derived from biological sources are free from infectious agents. When handling biological material appropriate safety precautions are recommended (Do not pipet by mouth; wear disposable gloves while handling biological material and performing the test; disinfect hands when finished the test). Biological material should be inactivated before disposal (e.g. in an autoclave). Disposables should be autoclaved or incinerated after use. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated areas swabbed with a suitable standard disinfectant or 70% alcohol. Material used to clean spills, including gloves, should be inactivated before disposal (e.g. in an autoclave).

Anti-HLA sera contain  $\text{NaN}_3$  as a preservative. The reagents contain < 0.1%  $\text{NaN}_3$  which is not considered to be a harmful concentration. Nevertheless avoid contact with the skin and mucous membranes. The copper and lead used in some plumbing systems can react with azides to form explosive salts. The quantities of azide used in this reagents are small; nevertheless when disposing of azide-containing materials, they should be flushed away with a large volume of water.

Disposal of all specimen and test materials should be in accordance with state and local law.

For quenching solution, Formaldehyde solution and Acridinorange / Ethidiumbromide (AO/EB) please adhere to the warnings and precautions of the manufacturer.

A yellow colouration of anti-HLA sera which still remains after thawing is indicating a change of the pH value. Those plates should **not** be used for the test.

Do not use **HISTO TRAY plates and rabbit complement** beyond the indicated expiration date on the label.

|                           |   |
|---------------------------|---|
| <b>Preservative:</b>      | < 0.1% $\text{NaN}_3$                             |
| <b>Storage:</b>           | ≤ - 20°C  |
| <b>Shelf life:</b>        | until the expiration date indicated on the labels |
| <b>Package / typings:</b> | according to information indicated on the kit     |

| Explanation of symbols used on Labelling  |                                 |   |  |
|---|---------------------------------|---|--|
|  | use by                          |  | Consult Instructions for use                         |
|  | Storage temperature             |  | sufficient for n tests                               |
| <b>ANTI-HLA-SERA</b>  | anti-HLA-Sera                   | <b>LOT</b>  | Batch code   |
| <b>COMPLEMENT RAB</b>   | rabbit complement               | <b>LYOPH</b>  | lyophilised  |
| <b>CONT</b>   | content, contains               | <b>MICROTESTTRAY</b>  | Microtest tray with predropped antisera and controls |
| <b>CONTROL -</b>  | negative Control                | <b>MONOCL</b>   | monoclonal   |
| <b>CONTROL +</b>  | positive Control                | <b>OR</b>   | or   |
| <b>HLA TYPING</b>   | Intended purpose:<br>HLA-typing | <b>POLYCL</b>   | polyclonal   |
| <b>HUM</b>  | Origin: human                   | <b>REF</b>  | Catalogue number                                     |
| <b>IFU</b>  | Instructions for use            | <b>WORKSHEET</b>  | Worksheet  |
| <b>IVD</b>  | For in vitro diagnostic use     |   |  |

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BAG Health Care GmbH