GEBRAUCHSINFORMATION

KANINCHENKOMPLEMENT

IN VITRO DIAGNOSTIKUM

REF 7018, 7023

Produktbeschreibung

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REF 7018, 7023

Description of product

The **rabbit complement** is produced of a large rabbit serum pool to ensure near identical complement activity in different lots.

Rabbit complement is e.g. used in the microlymphocytotoxicity test to detect HLA-antigens on the lymphocyte membrane.

Application

Dissolve lyophilized complement with 1 ml or 5 ml aqua dest. (according to package sizes). The reconstitution needs 10-15 minutes.

Reconstituted complement must be stored cool (2...8°C) and used within 3 - 4 hours.

DO NOT FREEZE dissolved rabbit complement!

Warnings and Precautions

The **rabbit complement** is designed for in vitro diagnostic use only and should be used by properly trained personnel.

Nevertheless all used biological material should be handled as potentially infectious, because no test method can guarantee that material derived from biological sources is free from infectious agents. When handling biological material appropriate safety precautions are recommended (do not pipet by mouth; wear disposable gloves while handling biological material and performing the test; disinfect hands when finished the test).

Biological material should be inactivated before disposal (e.g. in an autoclave). Disposables should be autoclaved or incinerated after use. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated areas swabbed with a suitable standard disinfectant or 70% alcohol. Material used to clean spills, including gloves, should be inactivated before disposal (e.g. in an autoclave). Disposal of all specimen and test materials should be in accordance with state and local law.

When used in accordance with the principles of Good Laboratory Practice, good standards of occupational hygiene and the instructions stated in this direction, the reagent supplied should not present a hazard to health.

Do not use **rabbit complement** beyond the indicated expiration date on the label.

Shelf life: lyophilized until the expiration date indicated on the label, dissolved see application

Storage: lyophilized 2...8°C or ≤ -20°C, dissolved see application

Packages: see information indicated on the label

Explanation of symbols used on Labelling			
	Use by		Consult Instructions for use
	Storage temperature		For in vitro diagnostic use
	Rabbit complement		Batch code
	Intended purpose: HLA-typing		Catalogue number
	lyophilised		or

Instructions for use in other languages see <http://www.bag-healthcare.com> or phone +49 (0)6404-925-125

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