Bemerkungen/Remarks:

| | | |
|---|---|-------------------------|
| wk = schwach/weak | + e = extra Reaktion/extra reaction | ↓ = carry over possible |
| [] = Split-Reaktion/Split-reactions | () = kann Reaktion zeigen/may show reactions | |
| * = Neues Antiserum/new antiserum Neue Position/new position | s.t.wk = manchmal schwach/sometimes weak | |
| | + = monoclonal | |

Bw4 und Bw6 Assoziationen Bw4 and Bw6 associations

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Directions for Use (NIH)

1. Bring the HISTO TRAY to a temperature of 18 - 22° C.
 2. Place 1 µl lymphocyte suspension (app. 2.000 - 3.000 cells) into each prealiquoted well.
 3. Incubate at a temperature of 18 - 22° C for 30 minutes.
 4. Add 5-6 µl rabbit complement.
 5. Incubate at a temperature 18 - 22° C for 60 minutes.*
 6. Add 3-4 µl Quenching solution.
 7. Read HISTO TRAY in a fluorescence microscope.
- * after point 5 you can proceed as follows:
6. Add 3-4 µl Eosin solution (5%) (soft touch method) and incubate for 5-10 minutes.
 7. Fix with 5-6 µl Formaldehyde solution (37%, pH 7.2) (soft touch method).
 8. Read after sedimentation of lymphocytes (at least 60 minutes).
 9. Cover the tray with a cover glass shortly before reading.

Rabbit Complement

Reconstitute lyophilized complement with 1 ml aqua dest.
Reconstituted complement use immediately. Do not freeze again.

Directions for Use Immuno beads (IMB) method

1. Isolate the lymphocytes by IMB technique.
2. Bring the HISTO TRAY to a temperature of 18 - 22° C.
3. Place 1 µl IMB-T-lymphocyte suspension (app. 1.000 cells) into each prealiquoted well.
4. Incubate at a temperature of 18 - 22° C for 30 minutes.
5. Add 5-6 µl rabbit complement (1.000 µl rabbit complement, 20 µl EB/AO).

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Microtest tray (72) with prealiquoted HLA-ABC anti sera and controls (human or monoclonal) for typing one specimen

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