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**Instructions for the Evaluation of
BAGene ABO-TYPE variant:**

1. Determination of mixes showing a positive reaction

Only those bands are classified as positive which have the correct size compared to a DNA length standard. The correct sizes of PCR products are given in the Worksheet as well as in the Evaluation Diagram (Table 1).

Example: positive band in mixes no. : 1, 2, 3, 4, 6, 8 and 11

**2. Checking positive reactions for ABO blood group genotyping using the
Evaluation Diagram (Table 2)**

Matching reaction pattern in the first column of Table 2. In Table 2 the reaction patterns are listed according to the lowest positive mix, so first you will find all reaction patterns that include mix 1, second all reaction patterns that include mix 2 and so on. Common genotypes are listed prior to rare variants.

Examples: 1+2+3+4+6+8+11 = O¹O²
 1+2+4+5+6+8+11 = O¹B¹
 1+2+4+6+7+8+11 = O¹A²
 1+4+6+8+11 = O¹O¹

**3. Confirmation and dDocumentation of reaction patterns on the worksheet of a
particular lot**

Read the results in the rows from left to right (red arrow on the worksheet). Fill in all positive and negative reactions writing + and – and compare the patterns with your results of the Evaluation Diagram (Table 2). An exception in the evaluation scheme are the mixes no. 10 and 11 denominated A^x and B³. These reactions do not indicate exclusively A^x and B³. They vary depending on which common alleles are present. Since the genotype can be indicated by positive reactions (band is present) and by negative reactions (band is missing), the differentiation of A^x und B³ has to be done with the complete reaction pattern.


Example: positive primer mixes
 1+2+3+4+6+8+11 = O¹O² (result acc. to evaluation diagram)
 complete reaction pattern = O¹O² (result acc. to worksheet)

Reaktionen-Nr. Reaction No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	Genotyp Genotype	Phänotyp Phenotype
	+	+	+	+	-	+	-	+	-	-	+	-	-	-	-	-	O ¹ O ²	O


4. Special remarks regarding genotyping the alleles A^x, A^w, B^x und B^w

Variant subgroup alleles are characterized by single nucleotide mutations (base substitutions). The same mutation can either indicate an A or a B subtype. Therefore the genotype is not only detected by a specific band in the corresponding primermix which indicates one particular mutation. In case of a positive „B-reaction“ (mix no. 5), the result is B^x oder B^w. In case of positivity of all „non-reactions“ (mixes no. 2, 4, 6 und 8), one and the same mutation leads to genotype A^x oder A^w. The reaction pattern of all PCR mixes must be considered completely.

Example:

A¹B^w

O¹A^w

5. Evaluation of results

The results of genotyping using **BAGene** ABO-TYPE variant have to be compared with serologic pre-testings, since the phenotype is decisive for the interpretation of variant subgroups.

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