Gebrauchsinformation
BAGene DNA-SSP Kits

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- ◆ **BAGene Partial D-TYPE** [20, 22, 24, 25, 28, 35-37], **Weak D-TYPE** [23, 25, 31], **D Zygoty-TYPE** [38], **KKD-TYPE** [56], **MNS-TYPE** [59], **HPA-TYPE** [66], **HNA-TYPE** [67]

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7. Literatur

• **BAGene ABO-TYPE**

1. **Ferguson-Smith MA**, Aitken DA, Turleau C, de Grouchy J: Localisation of the human ABO: Np-1: AK-1 linkage group by regional assignment of AK-1 to 9q34. *Hum Genet* 1976;34:35-43.
2. **Yamamoto F**, Clausen H, White T, Marken J, Hakamori S: Molecular genetic basis of the histo-blood group ABO system. *Nature* 1990; 345:229-33.
3. **Yamamoto F**, Hakamori S: Sugar-nucleotide donor specificity of histo-blood group A and B transferases is based on amino acid substitutions. *J Biol Chem* 1990;265: 19257-62.
4. **Yamamoto F**, McNeill PD, Yamamoto M, Hakomori S, Harris T, Judd WJ, Davenport RD: Molecular genetic analysis of the ABO blood group system: 1. Weak subgroups: A³ and B³ alleles. *Vox Sang* 1993;64:116-9.
5. **Yamamoto F**, McNeill PD, Yamamoto M, Hakomori S, Harris T: Molecular genetic analysis of the ABO blood group system: 3. A^x and B^(A) alleles. *Vox Sang* 1993;64:171-4.
6. **Chester MA**, Olsson ML. The ABO blood group gene: A locus of considerable genetic diversity. *Transfus Med Rev* 2001;15:177-200.
7. **Yip SP**. Sequence variation at the human ABO locus. *Ann Hum Genet* 2002;66:1-27.
8. **Gassner C**, Schmarda A, Nussbaumer W, Schönitzer D: ABO Glycosyltransferase Genotyping by polymerase chain reaction using sequence-specific primers. *Blood* 1996;88:1852-6.
9. **Müller TH**, Hallensleben M, Schunter F, Blasczyk R: Molekulargenetische Blutgruppendiagnostik. *Dt Ärztebl* 2001;98:A 317-322 [Heft 6].
10. **Seltsam A**, Hallensleben A, Kollmann A, Burkhart J, Blasczyk R. Systematic analysis of the ABO gene diversity within exons 6 and 7 by PCR-screening revealed new ABO alleles. *Transfusion* 2003;43:428-39.
11. **Olsson ML**, Irshaid NM, Hosseini-Maaf B, et al. Genomic analysis of clinical samples with serologic ABO blood grouping discrepancies: Identification of 15 novel A and B subgroup alleles. *Blood* 2001;98:1585-93.
12. **Seltsam A**, Hallensleben M, Kollmann A, Blasczyk R: The nature of diversity and diversification at the ABO locus. *Blood* 2003;102:3035-42.
13. **Seltsam A**, Das Gupta C, Wagner FF, Blasczyk R. Non-deletional ABO*O alleles express weak blood group A phenotypes. *Transfusion* 2005;45:359-65.

- **BAGene RH-TYPE, Partial D-TYPE, Weak D-TYPE, D Zygosity-TYPE**

14. **Chérif-Zahar B**, Localization of the human Rh blood group gene structure to chromosome region 1p34.3 – 1p36.1 by in situ hybridisation. *Hum Genet* 1991;86: 398-400.
15. **Blunt T**, Daniels G, Carritt B: Serotype switching in a partially deleted *RHD* gene. *Vox Sang* 1994;67:397-401.
16. **Rouillac C**, Colin Y, Hughes-Jones NC, et al.: Transcript analysis of D category phenotypes predicts hybrid Rh D-CE-D proteins associated with alteration of D epitopes. *Blood* 1995;85:2937-44.
17. **Rouillac C**, Le Van Kim C, Beolet M, Cartron J-P, Colin Y: Leu110Pro substitution in the RhD polypeptide is responsible for the DVII category blood group phenotype. *Am J Hematol* 1995;49:87-8.
18. **Mouro I**, Colin Y, Sistonen P, Le Pennec PY, Cartron J-P, Le Van Kim C: Molecular basis of the RhC^w (Rh8) and RhC^x (Rh9) blood group specificities. *Blood* 1995;86:1196-201.
19. **Beckers EAM**, Faas BHW, Simsek S, et al.: The genetic basis of a new partial D antigen: D^{DBT}. *Br J Haematol* 1996;93:720-7.
20. **Gassner C**, et al.: RHD/CE typing by polymerase chain reaction using sequence-specific primers *Transfusion* 1997;37:1020-6.
21. **Okada H**, Kawano M, Iwamoto S, Tanaka M, Seno T, Okubo Y, Kajii E: The *RHD* gene is highly detectable in RhD-negative Japanese donors. *J Clin Invest* 1997;100:373-379
22. **Wagner FF**, Gassner C, Müller TH, Schönitzer D, Schunter F, Flegel WA: Three molecular structures cause Rhesus D category VI phenotypes with distinct immunohematologic features. *Blood* 1998;91:2157-68.
23. **Legler TJ**, Maas JH, Blaschke M, Malekan H, Ohto R, Lynen R, Bustami N, Schwartz DWM, Mayr WR, Köhler M, Panzer S: *RHD* genotyping in weak D phenotypes by multiple polymerase chain reactions. *Transfusion* 1998;38:434-40.
24. **Omi T**, Takahashi J, Tsudo N, et al.: The genomic organization of the partial D category DVa : the presence of a new partial D associated with the DVa phenotype. *Biochem Biophys Res Commun* 1999;254:786-94.
25. **Wagner FF**, Gassner C, Müller TH, Schönitzer D, Schunter F, Flegel WA: Molecular basis of weak D phenotypes. *Blood* 1999;93:385-93.
26. **Flegel WA**, Wagner FF: Molecular genetics of *RH*. *Vox Sang* 2000;78(suppl 2): 109-15.
27. **Wagner FF**, Flegel WA: *RHD* gene deletion occurred in the *Rhesus box*. *Blood* 2000;95:3662-8.
28. **Avent ND**, Reid ME: The Rh blood group system: a review. *Blood* 2000;95:375-87.

29. **Legler TJ**, Wiemann V, Ohto H, Matuda I, Obara T, Uchikawa M, Köhler M: D^{Va} Category Phenotype and Genotype in Japanese Families. *Vox Sang* 2000;78:194-7.
30. **Wagner T**, Resch B, Legler TJ, Mossier C, Helmberg W, Köhler, M, Lanzer G: Severe HDN due to anti-Ce that required exchange transfusion. *Transfusion* 2000;40:571-4.
31. **Wagner FF**, Frohmajer A, Ladewig B, Eicher Ni, Lonicer CB, Muller TH, Siegel MH, Flegel WA: Weak D alleles express distinct phenotypes. *Blood* 2000;95:2699-708.
32. **Singleton BK**, Green CA, Avent ND, Martin PG, Smart E, Daka A, Narter-Olaga EG, Hawthorne LM, Daniels G: The presence of an *RHD* pseudogene containing a 37 base pair duplication and a nonsense mutation in Africans with the RhD-negative blood group phenotype. *Blood* 2000;95:12-18.
33. **Wagner FF**, Fromajer A, Flegel WA: *RHD* positive haplotypes in D negative Europeans. *BMC Genetics* 2001;2:10.
34. **Müller TH**, Hallensleben M, Schunter F, Blasczyk R: Molekulargenetische Blutgruppendiagnostik. *Dt Arztebl* 2001;98:A 317-322 [Heft 6].
35. **Shao CP**, Maas JH, Su YQ, Köhler M, Legler TJ: Molecular background of RhD-positive, D-negative, Del and weak D phenotypes in Chinese. *Vox Sang* 2002;83:156-61.
36. **Wagner FF**, Eicher NI, Jørgensen JR, Lonicer CB, Flegel WA: DNB: a partial D with anti-D frequent in Central Europe. *Blood* 2002;100:2253-6.
37. **Wagner FF**, Ladewig B, Angert KS, Heymann GA, Eicher NI, Flegel WA: The *DAU* allele cluster of the *RHD* gene. *Blood* 2002;100:306-11.
38. **Perco P**, Shao CP, Mayr WR, Panzer S, Legler TJ: Testing for the D zygosity with three different methods revealed altered Rhesus boxes and a new weak D type. *Transfusion* 2003;43:335-9.

- **BAGene KKD-TYPE**

39. **Lee S**, Naime DS, Reid ME, Redman CM: Molecular basis for the high-incidence antigens of the Kell blood group system. *Transfusion* 1997;37:1117-22.
40. **Lee S**, Wu X, Reid M, Zelinski T, Redman C: Molecular basis of the Kell (K1) phenotype. *Blood* 1995;85:912-6.
41. **Olives B**, Neau P, Bailly P, Hediger MA, Rousselet G, Cartron JP, Ripoche P: Cloning and functional expression of a urea transporter from human bone marrow cells. *J Biol Chem* 1994;269:31649.
42. **Lucien N**, Sidoux-Walter F, Olives B, Moulds J, Le Pennec PY, Cartron JP, Bailly P: Characterization of the gene encoding the human Kidd blood group/urea transporter protein. Evidence for splice site mutations in Jknull individuals. *J Biol Chem* 1998;273:12973-80.

43. **Cartron JP**, Ripoche P: Urea transport and Kidd blood groups. *Transfus Clin Biol* 1995;2:309.
44. **Olives B**, Merriman M, Bailly P, Bain S, Barnett A, Todd J, Cartron JP, Merriman T: The molecular basis of the Kidd blood group polymorphism and its lack of association with type 1 diabetes susceptibility. *Hum Mol Genet* 1997;6:1017.
45. **Neote K**, Mak JY, Kolakowski LFJ, Schall TJ: Functional and biochemical analysis of the cloned Duffy antigen: identity with the red blood cell chemokine receptor. *Blood* 1994;84:44-52.
46. **Chitnis CE**, Chaudhuri A, Horuk R, Pogo AO, Miller LH: The domain on the Duffy blood group antigen for binding *Plasmodium vivax* and *P. knowlesi* malarial parasites to erythrocytes. *J Exp Med* 1996;184:1531-6.
47. **Chaudhuri A**, Zbrzezna V, Polyakova J, Pogo AO, Hesselgesser J, Horuk R: Expression of the Duffy antigen in K562 cells. Evidence that it is the human erythrocyte chemokine receptor. *J Biol Chem* 1994;269:7835.
48. **Chaudhuri A**, Polyakova J, Zbrzezna V, Williams K, Gulati S, Pogo AO: Cloning of glycoprotein D cDNA, which encodes the major subunit of the Duffy blood group system and the receptor for the *Plasmodium vivax* malaria parasite. *Proc Natl Acad Sci USA* 1993;90:10793.
49. **Tournamille C**, Le van Kim C, Gane P, Cartron JP, Colin Y: Molecular basis and PCR-DNA typing of the Fya/fyb blood group polymorphism. *Hum Genet* 1995;95:407-10.
50. **Iwamoto S**, Omi T, Kajii E, Ikemoto S: Genomic organization of the glycoprotein D gene: Duffy blood group Fya/Fyb alloantigen system is associated with a polymorphism at the 44- amino acid residue. *Blood* 1995;85:622.
51. **Chaudhuri A**, Polyakova J, Zbrzezna V, Pogo AO: The coding sequence of Duffy blood group gene in humans and simians: restriction fragment length polymorphism, antibody and malarial parasite specificities, and expression in nonerythroid tissues in Duffy- negative individuals. *Blood* 1995;85:615.
52. **Mallinson G**, Soo KS, Schall TJ, Pisacka M, Anstee DJ: Mutations in the erythrocyte chemokine receptor (Duffy) gene: the molecular basis of the Fya/Fyb antigens and identification of a deletion in the Duffy gene of an apparently healthy individual with the Fy(a-b-) phenotype. *Br J Haematol* 1995;90:823-9.
53. **Tournamille C**, Colin Y, Cartron JP, Le van Kim C: Disruption of a GATA motif in the Duffy gene promoter abolishes erythroid gene expression in Duffy-negative individuals. *Nat Genet* 1995;10:224-8.
54. **Tournamille C**, Le van Kim C, Gane P, Le Pennec PY, Roubinet F, Babinet J, Cartron JP, Colin Y: Arg89Cys substitution results in very low membrane expression of the Duffy antigen/receptor for chemokines in Fy(x) individuals. *Blood* 1998;92:2147-56.
55. **Miller SA**, Dykes DD, Polesky HF: A simple salting out procedure for extracting DNA from human nucleated cells. *Nucleic Acids Res* 1988;16:1215.

56. **Rožman P**, Dovč T, Gassner C: Differentiation of autologous ABO, RHD, RHCE, KEL, JK, and FY blood group genotypes by analysis of peripheral blood samples of patients who have recently received multiple transfusions. *Transfusion* 2000;40:936-42.

• **BAGene MNS-TYPE**

57. **Landsteiner K**, Levine P: A new agglutinable factor differentiating individual human bloods. *Proc Soc exp Biol* 1927;24:600.)

58. **Walsh RJ**, Montgomery C: A new human isoagglutinin subdividing the MN blood groups. *Nature* 1947;160:504.

59. **Shih MC**, Yang LH, Wang NM, Chang JG: genomic typing of human red cell Miltenberger glycoporphins in a Taiwanese population. *Transfusion* 2000;40:54-61.

60. **Storry JR**, Reid ME, Fetcs S, Huang CH: Mutations in *GYPB* exon 5 drive the S-s-U+^{var} phenotype in persons of African descent: implications for transfusion. *Transfusion* 2003;43:1738-47.

• **BAGene HPA-TYPE**

61. **Ballem PJ**, Buskard NA, Decary F, et al.: Post-transfusion purpura secondary to passive transfer of anti-PI^{A1} by blood transfusion. *Br J Hematol* 1987;66:113-4.

62. **Mueller-Eckhardt C**, Kiefel V, Grubert A, et al.: 348 cases of suspected neonatal alloimmune thrombocytopenia. *Lancet* 1989;1:363-6.

63. **Panzer S**, Kiefel V, Bartram CR, et al.: Immune thrombocytopenia more than a year after allogeneic marrow transplantation due to antibodies against donor platelets with anti-PI^{A1} specificity: evidence for a host derived immune reaction. *Br J Hematol* 1989;71:259-64.

64. **Mueller-Eckhardt C**, Kiefel V, Santoso S: Review and update of platelet alloantigen systems. *Transfus Med Rev* 1990;4:98-109.

65. **Santoso S**, Kiefel V: Human platelet specific alloantigens: update. *Vox Sang* 1998;74 (Suppl):249-53.

66. **Jau-Yi L**, Ying-Ju C, Hui-Yu H, Jeong-Shi L, Cheng-Hwai T: PCR with sequence-specific primer-based simultaneous genotyping of human platelet antigen-1 to -13w. *Transfusion* 2002;42:1089-95.

• **BAGene HNA-TYPE**


67. **Carl B**, Kroll H, Bux J, Bein G, Santoso S: B-lymphoblastoid cell lines as a source of reference DNA for human platelet and neutrophil antigen genotyping. *Transfusion* 2000;40:62-68.

68. **Sachs UJH**, Chavakis T, Fung L, Lohrenz A, Bux J, Reil A, Ruf A, Santoso S: Human alloantibody anti-Mart interferes with Mac-1-dependent leukocyte adhesion. *Blood* 2004;104,3:727-34.
69. **Sachs UJH**, Reil A, Bauer C, Bux J, Bein G, Santoso S: Genotyping of human neutrophil antigen-5a (Ond). *Transfus Med* 2005;15:115-117.

- **Allgemein**

70. **Maniatis et al.**, 1989. *Molecular Cloning: A Laboratory Manual*.
New York: Cold Spring Harbour Laboratory
71. **Beutler E et al.**, 1990. *BioTechniques* 9:166
72. **Olerup O**, Zetterquist H. HLA-DR typing by PCR amplification with sequence-specific primers (PCR-SSP) in 2 hours: an alternative to serological DR typing in clinical practice including donor-recipient matching in cadaveric transplantation. *Tissue Antigens* 1992;39:225-35.
73. **Olerup O**, Zetterquist H. DR "low-resolution" PCR-SSP typing - a correlation and an update. *Tissue Antigens* 1993;1:55-6.
74. **Chen EY**, Liao YC, Smith DH, Barrera-Saldana HA, Gelinis RE, Seeburg PH. The human growth hormone locus: nucleotide sequence, biology, and evolution. *Genomics* 1989;4:479.

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BAG Health Care GmbH