## NEU: BAG-ProSpore Dry Heat

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10/09

## PERFORMANCE DATA for ProSpore-DH<sup>®</sup>

**ORGANISM:** *Bacillus atrophaeus* derived from (ATCC #9372)

**POPULATION:\***  $1.9 \times 10^6$  CFU'S\*\* / ampoule

**RESISTANCE CHARACTERISTICS:**

**D<sub>160</sub> VALUE:\*\*\*** **4.2** minutes (Dry Heat, 160°C – This accuracy shall not exceed +/- 0.2)

\*After preliminary heat treatment

\*\* Colony forming units

\*\*\* Determined on Primary Spore Crop using Fraction-Negative after graded exposures

AGENT	CONDITIONS	SURVIVES	KILLED
Dry Heat	160 ± 2°C	18.0 min	43.2 min

**INCUBATION:** 4 days at 37°C +/- 2°C in supplied Tryptic Soy Broth with Phenol Red  
**STORAGE:** 15-27°C (60-80°F), 30-70% RH, away from sterilizing agents, direct sunlight and all other forms of UV light. (DO NOT REFRIGERATE)  
**DISPOSAL:** Do not use after expiration date. Sterilize all cultures before discarding.

**05/10**

**D100**

### ProSpore DH - Instructions for Use

1. Include several biological indicators with each sterilizer cycle to be monitored. Place in the most difficult to sterilize location in chamber
2. After sterilization:
  - a) Let the ampoule cool for at least ten (10) minutes
  - b) Prior to placing in the workstation (laminar flow hood) disinfect the ampoule with 70% isopropanol alcohol
  - c) Snap the top off the ampoule by using the thumb as a hinge and applying pressure at the scored glass line
  - d) Slowly and carefully pour the contents of the ampoule into a test tube of Modified Tryptic Soy Broth, making sure none of the sand is stuck to the side of the test tube. One tube of prepared media is provided with each ampoule.
  - e) Incubate at 37 +/- 2°C
  - f) Observe tubes daily for growth
  - g) Final results are read at 96 hours
3. Positive control:
  - a) Mark one (1) tube to be positive control.
  - b) Repeat steps b) through f) from Section 2 (above)
4. Record Results. Results are valid only if control ampoule shows growth.