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*Microtest trays (60) with prealiquoted antisera and controls for HLA-B27-typing*

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### Directions for Use (NIH)

1. Bring the HISTO TRAY to a temperature of 18 - 22° C.
2. Place 1 µl lymphocyte suspension (app. 2.000 - 3.000 cells) into each prealiquoted well.
3. Incubate at a temperature of 18 - 22° C for 30 minutes.
4. Add 5-6 µl rabbit complement.
5. Incubate at a temperature 18 - 22° C for 60 minutes.\*
6. Add 3-4 µl Quenching solution.
7. Read HISTO TRAY in a fluorescence microscope.

\* after point 5 you can proceed as follows:

6. Add 3-4 µl Eosin solution (5%) (soft touch method) and incubate for 5-10 minutes.
7. Fix with 5-6 µl Formaldehyde solution (37%, pH 7.2) (soft touch method).
8. Read after sedimentation of lymphocytes (at least 60 minutes).
9. Cover the tray with a cover glass shortly before reading.

### Rabbit Complement

Reconstitute lyophilized complement with 1 ml aqua dest. Reconstituted complement use immediately. Do not freeze again.

### Directions for Use Immuno beads (IMB) method

1. Isolate the lymphocytes by IMB technique.
2. Bring the HISTO TRAY to a temperature of 18 - 22° C.
3. Place 1 µl IMB-T-lymphocyte suspension (app. 1.000 cells) into each prealiquoted well.
4. Incubate at a temperature of 18 - 22° C for 30 minutes.
5. Add 5-6 µl rabbit complement (1.000 µl rabbit complement, 20 µl EB/AO).
6. Incubate for 60 minutes at a temperature of 18 - 22° C in darkness.
7. Add 3-4 µl diluted Quenching-solution (2.000 µl Quenching + 1.000 µl EDTA 8%).
8. Read HISTO TRAY in a fluorescence microscope.