**DE**

## Gebrauchsinformation

### TERASAKI Park MEDIUM

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## Instructions for Use

### TERASAKI Park MEDIUM

Cell culture medium for storage of cells (e.g. lymphocytes)

**REF 7028**

#### Description of product

The intended use of **TERASAKI Park MEDIUM** is the storage of cells, e.g. lymphocytes. The solution is ready to use. It is suitable to maintain physiological conditions for cells, e.g. human lymphocytes.

In **TERASAKI Park MEDIUM**, isolated cells (e.g. lymphocytes) may be kept up to 4 days without noticeable loss of viability.

#### Application e.g. for lymphocytes

Isolate lymphocytes as usual. After the last washing step, the sediment is topped up with 5-10 ml of **TERASAKI Park MEDIUM** and stored refrigerated at 2...8°C.

In order to receive the required amount of cells for e.g. HLA-typing proceed as follows:

1. Centrifuge the cell suspension for 10 minutes at 110 x g.
2. Discard the supernatant.
3. Adjust the sediment with **TERASAKI Park MEDIUM** to the required amount of cells (e.g. 2000-3000 cells/µl).

#### Warnings and Precautions

Do not use **TERASAKI Park MEDIUM** in case of changing the colour from orange to yellow. Do not use **TERASAKI Park MEDIUM** beyond the indicated expiration date on the label.

**Shelf life:** until the expiration date indicated on the label, after opening to be used within 1-2 weeks

**Storage:** 2...8°C

**Package:** 100 ml

Explanation of symbols used on Labelling	
	Storage temperature
	Use by
	Consult Instructions for use
<b>CELLCULTURE MEDIUM</b>	Cellculture medium
<b>CONT FCS 0.5%</b>	Contains 0.5% foetal calf serum
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number

**Version: 2/2015 / Issue: 2015-08**