DE

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Information: Instructions for Use

Mastermixes / Magnesium Chloride for HISTO SPOT[®] SSO Kits

Test kits for tissue typing of HLA alleles on a molecular genetic basis

IVD

CE 0123

- REF 726080: HISTO SPOT[®] Mastermix A** (for HISTO SPOT[®] A 4D, **REF 726010**)
REF 726081: HISTO SPOT[®] Mastermix B (for HISTO SPOT[®] B 4D, **REF 726020**)
REF 726083: HISTO SPOT[®] Mastermix C (for HISTO SPOT[®] C 4D, **REF 726030**)
REF 726084: HISTO SPOT[®] Mastermix DRB1 (for HISTO SPOT[®] DRB1 4D, **REF 726040**)
REF 726204: HISTO SPOT[®] Mastermix DRB3/4/5 (for HISTO SPOT[®] DRB3/4/5, **REF 726045**)
REF 726082: HISTO SPOT[®] Mastermix DQ (for HISTO SPOT[®] DQB1 4D/ DQA1, **REF 726051**)
REF 726086: HISTO SPOT[®] Mastermix DPB1 (for HISTO SPOT[®] DPB1, **REF 726060**, **726061**)
REF 726203: HISTO SPOT[®] Mastermix Coeliac Disease (for HISTO SPOT[®] Coeliac Disease, **REF 726071**)
REF 726087: HISTO SPOT[®] Magnesium Chloride (for all HISTO SPOT Typing Kits)

Mastermix I X 1100 µl **Mastermix**, ready to use, contains biotinylated primers for the chosen locus, dNTPs, Taq polymerase, reaction buffer, 0.05% sodium azide

MgCl₂ **Magnesium Chloride**, 6 mM, ready to use, contains 0.001 % Proclin[®] 300

X = A, B, C, DRB1, **DRB3/4/5**, DQB1, DQ, **DPB1** or **Coeliac**

The Mastermixes and the magnesium chloride are only suitable for use with HISTO SPOT[®] SSO-Kits. The Instructions for Use for HISTO SPOT[®] SSO-Kits must be observed. The Instructions for Use is in each HISTO SPOT[®] SSO kit on CD and can be downloaded from <http://www.bag-healthcare.com/en/Diagnostika/Downloads/>

STORAGE

The mastermixes and the magnesium chloride should be stored at 2...8°C. The expiry date is indicated on the label and is valid for the originally sealed reagent. The opened reagent should be used within 3 months.

WARNINGS AND PRECAUTIONS

HISTO SPOT[®] is designed for in vitro diagnostic use and should be used by properly trained, qualified staff. The magnesium chloride solution contains ProClin[®]300 in a concentration of 0.001% as preservative. The mastermixes contain the preservative sodium azide. The reagents contain < 0.1% sodium azide which is not considered to be a harmful concentration. Nevertheless avoid contact with the skin and mucous membranes. Sodium azide may react with lead and copper plumbing to form explosive metal azides. While disposing of sodium azide containing solutions down laboratory sinks, flush the drains with a large volume of water to prevent azide build-up.

All work with reagents should be handled with the appropriate precautions. Wear eye protection, laboratory coats and disposable gloves when handling the reagents. Avoid contact of these materials with the skin, eyes or mucous membranes. If contact does occur, immediately wash with large amounts of water.

If spills of reagents occur, dilute with water before wiping dry.

Disposal should be in accordance with country, federal, state and local regulations.

Material Safety Data Sheets (MSDS) are available to download at www.bag-healthcare.com.

EXPLANATION OF SYMBOLS USED ON LABELING

	For in vitro diagnostic use		Use by		Batch code
	Storage temperature		Catalogue number		Consult instructions for use
Mastermix A	Mastermix for amplification of the HLA-A locus				
Mastermix B	Mastermix for amplification of the HLA-B locus				
Mastermix C	Mastermix for amplification of the HLA-C locus				
Mastermix DRB1	Mastermix for amplification of the HLA-DRB1 locus				
Mastermix DRB3/4/5	Mastermix for amplification of the HLA-DRB3/4/5 loci				
Mastermix DQ	Mastermix for amplification of the HLA-DQB1 and DQA1 locus				
Mastermix DPB1	Mastermix for amplification of the HLA-DPB1 locus				
Mastermix Coeliac	Mastermix for amplification of specific HLA-DQ loci				
MgCl₂	Magnesium chloride solution				

Version: 6/2014

Changes to version 5/2014 are highlighted

Issued: 2014-10