## GEBRAUCHSINFORMATION

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FOR IN VITRO DIAGNOSTIC USE

#### **1. Product description**

Anti-c and Anti-E monoclonal (IgM) are prepared from monoclonal, human IgM antibodies. The clone numbers are given on the labels of the test reagents.

Anti-c and Anti-E monoclonal (IgM) are designed for use in tube and plate tests, and provides a specific, qualitative test for the detection of the corresponding antigens on human red blood cells.

The test reagents contain < 0.1% NaN<sub>3</sub> as preservative.

#### **2. Biological principle of the test**

The test used with this blood grouping reagents is based on the principle of hemagglutination. Incubation of test red cells with the monoclonal test reagents will result in a specific antigen-antibody reaction if the corresponding antigen is present on the test cells. Visible detection of this reaction is demonstrated by agglutination of the cells. No agglutination indicates a negative test result, and within the accepted limitations of the test procedure, indicates the absence of the corresponding antigen.

#### **3. Storage and Shelf Life**

Store the test reagents at 2...8°C. Do not freeze! Allow the test reagents to reach room temperature (18...25°C) before use. Return reagents to 2...8°C for storage as appropriate, immediately after use. After opening the bottle the test reagents can be used until the expiry date printed on the label, if appropriate storage conditions be observed. Do not use the reagents after the expiry date printed on the label.

#### **4. Specimen preparation**

Blood samples should be collected by approved medical procedure. Blood collected without or with anticoagulant (EDTA, heparin) is acceptable. Do not use haemolytic samples. Testing should be performed without delay if possible. If this is not possible, store blood samples at 2...8°C.

Prolonged storage of red cells prior to testing may result in deterioration of red cell antigens and resultant weaker than expected test reactions (s. 9. Important Directions / Limitations of Procedure).

#### **5. Additional Materials Required**

Isotonic saline, test plates for Blood Group Typing, test tubes (75 x 12 mm), disposable Pasteur Pipettes, centrifuge

#### **6. Test procedure**

##### **Plate test**

1. Prepare a 10% suspension of test red cells in isotonic saline.
2. Place 1 drop of test reagent and 1 drop of the prepared suspension of test red cells on a test plate and mix.
3. Incubate for 5 - 10 minutes at room temperature.
4. By slowly rotation of the plate, examine macroscopically for agglutination.

##### **Tube test**

1. Wash the test red cells at least once with isotonic saline and prepare an approximately 3% suspension of test red cells in isotonic saline.
2. Place 1 drop of test reagent and 1 drop of the prepared suspension of test red cells into a labelled test tube, mix and incubate for 10 - 15 minutes at room temperature.
3. Centrifuge 1 minute at 400 x g (1500 rpm) or at an alternative rpm with an appropriate time adjustment.
4. Resuspend the cells by gently shaking the tube and examine macroscopically for agglutination.

##### **Directions:**

Do not examine tests microscopically.

Red blood cell suspensions known to be positive (ideally heterozygous cells) and negative for the antigen, Rh-control for monoclonal test reagents and a patient control should always be included in the test.

Use at least two different test reagents to determine the antigen. By use of two monoclonal test reagents two different clones should be used.

#### **7. Interpretation of test results**

Agglutination of test red cells with the test reagent indicates the presence of the corresponding antigen (within the accepted limitations of the test procedure).

No agglutination of test red cells with the test reagent indicates the absence of the corresponding antigen (within the accepted limitations of the test procedure).

If no agglutination occurs with the test red cells known to be positive for the antigen or if agglutination occurs with the test red cells known to be negative for the antigen or the Rh control for monoclonal test reagents or with the patient control the test results should not be interpreted.

If different test results occur with two different test reagents, repeat the determination of the antigen with another test method and/or another test reagent.

Pay attention to the limitations of procedure and important directions.

(s. 9. Important Directions / Limitations of Procedure).

### **8. Stability Of The Reaction**

All test results should be interpreted immediately upon completion of the test.

### **9. Important Directions / Limitations of Procedure**

1. The test reagents are designed for in vitro diagnostic use only and should be used by properly trained individuals.
2. On rare occasion, red cells coated in vivo with immunoglobulin may agglutinate spontaneously and non-specifically. In such instances similar phenomena would most likely occur in the ABO grouping test and blood grouping tests of other blood group systems as well. Rh control for monoclonal test reagents and patient autologous serum are suitable controls. If the control test yields a positive reaction, a valid interpretation of the blood typing result cannot be made.
3. The use of unwashed test red cells suspended in plasma or serum may promote false positive reactions such as those associated with rouleaux formation, or autoantibodies. The use of well washed red cells may reduce the incidence of such false positive reactions.
4. The strength of positive reactions depends on the age of the used blood.
5. Light cloudiness does not influence the reactivity of the product.
6. If the plate test is read too late, the appearance caused by drying may simulate false positive results. Therefore the plate test results should be interpreted after 10 minutes, at most 15 minutes.
7. False negative results or unexpected weak reactions may be caused by an insufficient cell concentration, insufficient incubation temperature or time and/or insufficient centrifugation, but also by storing the red cells for too long and/or under inappropriate conditions.
8. Other variables such as improper technique, inappropriate centrifugation or incubation, improperly cleaned glassware, incorrect saline pH and/or contaminated materials and samples may cause false negative or false positive results.
9. Microbiological contamination of the test reagents must be avoided as this may reduce the life of the products and cause erroneous results.
10. Some red cells may express quantitatively weak or Rh variants and may therefore demonstrate weaker than expected reactions with the test reagents, in the plate test also no reactions. Further clarification and specification of the result can be carried out with **BAGene** (SSP-Kits for the determination of Rh attributes on a molecular genetic basis).
11. No single centrifugation speed or time can be recommended for all types of available centrifuges or test applications. Centrifuges should be calibrated individually to determine the optimal time and speed required to produce a clear supernatant and a clearly delineated red cell button that can be easily resuspended.
12. Deviation from the recommended Instructions for Use may result in less than optimal product performance. User-defined deviations such as modifications of test procedures or freezing of serum on microtiter plates etc. may require validation by the user.
13. For interpretation of the test results, consider if transfusion or transplantation had happened. Take the case history of the transfusion or transplantation and also the patient's medication history into consideration.

### **10. Warnings and Precautions**

Human source material used to produce these reagents has been tested and found negative for HBsAg and HIV and HCV antibodies. Nevertheless all used biological material, especially the test red cells, should be handled as potentially infectious, because no test method can guarantee that material derived from biological sources are free from infectious agents. When handling biological material appropriate safety precautions are recommended (Do not pipette by mouth; wear disposable gloves while handling biological material and performing the test; disinfect hands when finished the test).

Biological material should be inactivated before disposal (e.g. in an autoclave). Disposables should be autoclaved or incinerated after use. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated areas swabbed with a suitable standard disinfectant or 70% alcohol. Material used to clean spills, including gloves, should be inactivated before disposal (e.g. in an autoclave).

The test reagents contain  $\text{NaN}_3$  as a preservative. The reagents contain  $< 0.1\%$   $\text{NaN}_3$  which is not considered to be a harmful concentration. Nevertheless avoid contact with the skin and mucous membranes. The copper and lead used in some plumbing systems can react with azides to form explosive salts. The quantities of azide used in this reagents are small; nevertheless when disposing of azide-containing materials, they should be flushed away with a large volume of water.



Disposal of all specimen and test materials should be in accordance with state and local law.

**11. Packages:** see price list

### **12. References**

Applied Blood Group Serology, PD Issitt and DJ Anstee, 4<sup>th</sup> Edition, Montgomery Scientific, Durham SC, 1998

Technical manual of the American Association of Blood Banks, 15<sup>th</sup> ed., 2005

Erklärung der Symbole auf den Etiketten / <i>Explanation of symbols used on Labelling</i>		
 <b>IVD</b>	In-vitro-Diagnostikum / <i>For in vitro diagnostic use</i>	 <b>MONOCL IGM</b> Monoklonal IgM / <i>Monoclonal IgM</i>
	Lagertemperatur / <i>Storage temperature</i>	 <b>CLONE</b> Klon / <i>Clone</i>
 <b>LOT</b>	Lot-Nr. / <i>Batch code</i>	 <b>ORIG HUM</b> Ursprung: human / <i>Origin: human</i>
	Verwendbar bis / <i>Use by</i>	 <b>CONT NaN<sub>3</sub></b> Enthält Natriumazid / <i>Contains Natriumazide</i>
 <b>REF</b>	Bestell-Nr. / <i>Catalogue number</i>	 <b>TIT</b> Titer / <i>Titer</i>
	Gebrauchsinformation beachten / <i>Consult instructions for use</i>	

**BAG Health Care GmbH**

Amtsgerichtsstraße 1-5  
35423 Lich/Germany

Tel.: +49 (0) 6404 / 925 - 0

Fax: +49 (0) 6404 / 925 - 250

[www.bag-healthcare.com](http://www.bag-healthcare.com)

[info@bag-healthcare.com](mailto:info@bag-healthcare.com)

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Tel.: +49 (0) 6404 / 925 - 450

Fax: +49 (0) 6404 / 925 - 460

[verkauf@bag-healthcare.com](mailto:verkauf@bag-healthcare.com)

**Customer Service:**

Tel.: +49 (0) 6404 / 925 - 125

Fax: +49 (0) 6404 / 925 - 421

[service@bag-healthcare.com](mailto:service@bag-healthcare.com)