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06/12

PERFORMANCE DATA for DriAMP™

ORGANISM: *Bacillus atrophaeus derived from (ATCC #9372*)*

POPULATION:* 3.4×10^6 CFU'S** / ampoule

RESISTANCE CHARACTERISTICS:

D₁₆₀ VALUE:*** **4.0** minutes

* Culture is traceable to a recognized culture collection identified in USP and ISO 11138

** Colony forming units

*** D-values are calculated using the Limited-Holcomb-Spearmen-Karber method.

AGENT	CONDITIONS	SURVIVES	KILLED
Dry Heat	160 ± 2°C	18.1 min	42.1 min

INCUBATION: 72 hours at 36°C - 38°C in supplied Releasat® Culture Medium

STORAGE: 15-27°C (60-80°F), 30-70% RH, away from sterilizing agents, direct sunlight and all other forms of UV light. (DO NOT REFRIGERATE)

DISPOSAL: Do not use after expiration date. Sterilize all cultures before discarding.

05/13

DH-105

DriAMP - Instructions for Use

1. Include several biological indicators with each sterilizer cycle to be monitored. Place in the most difficult to sterilize location in chamber.
2. After sterilization:
 - a. Let the ampoule cool for at least ten (10) minutes before culturing.
 - b. Before the ampoule is placed in the workstation, the ampoule should be sprayed with 70% IPA.
 - c. Remove the caps of the ampoules as follows:
 - Hold the top of the ampoule with one hand and the body of the ampoule with the other hand.
 - Place thumbs tip-to-tip over the scored glass line at the neck of the ampoule.
 - Snap the top off of the ampoule by using the thumb tips as a hinge and applying pressure at the scored glass line.
 - d. Slowly pour the contents from the ampoule into a test tube of Releasat® medium, one tube for each ampoule.
 - e. Incubate at 36-38°C.
 - f. Results are read at 72 hours.
 - g. Observe tubes daily for growth.
Yellow = growth = non-sterile / Red = no growth = sterile.
3. **Positive Control:** One or more positive controls should be included in each test series. Transfer a non-sterilized ampoule to Releasat® culture medium and incubate with test series. Yellow color and turbidity indicates that the medium possesses suitable growth promoting qualities and that the carrier contains viable spores. If positive control does not grow, do not use the units from that package.