

BAG - Bowie & Dick Test Systems



Lead - free

BAG-BD-Check II

Disposable Bowie-Dick test pack.

Pre-assembled test pack for the daily steam penetration test in prevacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7563

Packaging unit: 30 tests/box



BAG-BD-Check

Disposable Bowie-Dick test pack.

Pre-assembled test pack for the daily steam penetration test in prevacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7523

Packaging unit: 30 tests/box



Lead - free

BAG-AutoCheck Kit II

Reusable Bowie-Dick test pack. Use with small CCS holder.

Paper test packs and chemical indicators for the daily steam penetration test in prevacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7561

Packaging unit: 30 tests/box



BAG-AutoCheck Kit

Reusable Bowie-Dick test pack. Use with small CCS holder.

Paper test packs and chemical indicators for the daily steam penetration test in prevacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7521

Packaging unit: 30 tests/box



BAG-GreenCard *(available on request only)*

Single-use Bowie-Dick test pack. Use with GreenCard holder. Test card the size of credit card for the daily steam penetration test in prevacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 9 x 5 cm, 1 mm thin

Lead - free

Order-No.: 7564

Packaging unit: 30 tests/box

BAG-BD-Kit

Reusable Bowie-Dick test pack. Use with large CCS holder. Paper test packs and chemical indicators for the daily steam penetration test in prevacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 28 x 22 cm

BAG-BD-Kit II (Lead-free)

Order-No.: 7566

Packaging unit: 60 tests/box

BAG-BD-Kit

Order-No.: 7526

Packaging unit: 60 tests/box

BAG-BD-Kit 121 °C

Order-No.: 7528

Packaging unit: 60 tests/box

BAG-AutoCheck II / BAG-AutoCheck

Indicator sheets for the daily steam penetration test acc. to ISO 11140-3. Size: 28 x 22 cm

BAG-AutoCheck II (Lead-free)

Order-No.: 7566

Packaging unit: 100 indicators/box

BAG-AutoCheck

Order-No.: 7526

Packaging unit: 100 indicators/box

BAG-CCS test holder

For use with test packs 14 x 11 cm (BAG-AutoCheck-Kits, BAG-Biocheck-Kit)

Order-No.: 74651

Packaging unit: 1

Bowie-Dick Test holder

For use with test packs 28 x 22 cm (BAG-BD-Kits)

Order-No.: 7460

Packaging unit: 1

BAG-GreenCard holder

For use with BAG-GreenCard

Order-No.: 7565

Packaging unit: 1

05/08