

BAGene

HPA-TYPE

CE

LOT 0806 HP

IVD

2009-12

REF

6660

Worksheet und Auswertetabelle / Worksheet and Evaluation diagram

Reaktions-Nr. / Reaction No.	1	2	3	4	5	6	7	8	9	10	11	12
PCR-Produkt (Größe in bp) PCR product (size in bp)	196	196	241	241	230	230	252	252	246	246	233	233
HPA-Merkmale HPA specificities	HPA 1a	HPA 1b	HPA 2a	HPA 2b	HPA 3a	HPA 3b	HPA 4a	HPA 4b	HPA 5a	HPA 5b	HPA 15a	HPA 15b
	Beispiele / Examples											
⇒ HPA - 1a/b; 2a/b; 3a/b; 4a/b; 5a/b; 15a/b	+	+	+	+	+	+	+	+	+	+	+	+
HPA - 1a/a; 2a/a; 3a/a; 4a/a; 5a/b; 15a/a	+	-	+	-	+	-	+	-	+	+	+	-
HPA - 1b/b; 2b/b; 3b/b; 4b/b; 5b/b; 15b/b	-	+	-	+	-	+	-	+	-	+	-	+
⇒ HPA - 1a/b; 2a/a; 3a/b; 4a/a; 5b/b; 15a/b	+	+	+	-	+	+	+	-	-	+	+	+

Ergebnis Result	1	2	3	4	5	6	7	8	9	10	11	12	Genotyp Genotype

Proben-ID / Sample-ID:

Name:

Geb.-Datum / Birthdate:

Ergebnis / Result:

Datum / Date:

Unterschrift / Signature:

Gelbild / Gel Picture

VER 06/2008