

BAG - Bowie & Dick Test Systems



Lead - free

BAG-BD-Check II

Disposable Bowie-Dick test pack.

Pre-assembled test pack for the daily steam penetration test in pre-vacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7563

Packaging unit: 30 tests/box



BAG-BD-Check

Disposable Bowie-Dick test pack.

Pre-assembled test pack for the daily steam penetration test in pre-vacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7523

Packaging unit: 30 tests/box



Lead - free

BAG-AutoCheck Kit II

Reusable Bowie-Dick test pack. Use with small CCS holder.

Paper test packs and chemical indicators for the daily steam penetration test in pre-vacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7561

Packaging unit: 30 tests/box



BAG-AutoCheck Kit

Reusable Bowie-Dick test pack. Use with small CCS holder.

Paper test packs and chemical indicators for the daily steam penetration test in pre-vacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4. Size: 14 x 11 cm

Order-No.: 7521

Packaging unit: 30 tests/box



Lead - free

BAG-GreenCard *(available on request only)*

Single-use Bowie-Dick test pack. Use with GreenCard holder.
 Test card the size of credit card for the daily steam penetration test in pre-vacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4.
 Size: 9 x 5 cm, 1 mm thin

Order-No.: 7564

Packaging unit: 30 tests/box

BAG-BD-Kit

Reusable Bowie-Dick test pack. Use with large CCS holder.
 Paper test packs and chemical indicators for the daily steam penetration test in pre-vacuum steam sterilizers used as an alternative to the standard towel pack acc. to ISO 11140-4.
 Size: 28 x 22 cm

BAG-BD-Kit II (Lead-free)

Order-No.: 7566

Packaging unit: 60 tests/box

BAG-BD-Kit

Order-No.: 7526

Packaging unit: 60 tests/box

BAG-BD-Kit 121 °C

Order-No.: 7528

Packaging unit: 60 tests/box

BAG-AutoCheck II / BAG-AutoCheck

Indicator sheets for the daily steam penetration test acc. to ISO 11140-3. Size: 28 x 22 cm

BAG-AutoCheck II (Lead-free)

Order-No.: 7566

Packaging unit: 100 indicators/box

BAG-AutoCheck

Order-No.: 7526

Packaging unit: 100 indicators/box

BAG-CCS test holder

For use with test packs 14 x 11 cm (BAG-AutoCheck-Kits, BAG-BioCheck-Kit)

Order-No.: 74651

Packaging unit: 1

Bowie-Dick Test holder

For use with test packs 28 x 22 cm (BAG-BD-Kits)

Order-No.: 7460

Packaging unit: 1

BAG-GreenCard holder

For use with BAG-GreenCard

Order-No.: 7565

Packaging unit: 1